## 1. Product Name

Super Safe Plastic Fibre Optics

## 2. Product Code

17903 - 1m x 100 strands  
19869 - 1m x 150 strands  
17905 - 2m x 100 strands  
17907 - 2m x 150 strands

17908 - 2m x 200 strands  
19870 - 3m x 150 strands  
19871 - 3m x 200 strands

## 3. Colour

Strands are transparent until the coloured light from the light source illuminates the strands

## 4. Brief Description

The ultimate in fibre optic safety!

Use of all fibre optic products should be supervised. Whilst it is NOT recommended that these are mouthed, plastic tails are the safest option available for those who may do so.

Requires light source 19210.

## **5. Contents**

- 1 fibre optic harness

## **6. Snoezelen® Stimulations**

- Sight
- Touch

## **7. Best Use**

Use with a ROMPA® light source in a darkened environment. For interactive use, use with one of ROMPA®'s Wi Fi controllers or with Sensory Magic™.

## **8. Compatible ROMPA® Products**

19210 Light Source by ROMPA®  
20772 Fibre Optic Wooden 'Comb'  
17164 Fibre Optic Comb  
19545 Curved Fibre Optic Comb

*Please refer to the separate instructions for the Light Source*

## **9. Starting Up**

- Remove the packaging and temporary ties.
- Attach the boss to the light source.

## **10. Detailed Description**

The ultimate in fibre optic safety!

Use of all fibre optic products should be supervised. Whilst it is NOT recommended that these are mouthed, plastic tails are the safest option available for those who may do so. Requires light source 19210. Robust, the strands are ultra-resistant to accidental damage and are equally bright as standard fibre optics. Durable and long-lasting

## **SUGGESTIONS FOR USE OF FIBRE OPTIC PRODUCTS**

Because of the light source, fibre optic sprays change colour constantly along their length. The changing colours are visually stimulating and encourage focusing, yet are calming and relaxing too. The strands provide fascinating tactile stimulation too. Some suggestions for use are:

1. Attach a light source to a wall with the spray cascading down like a waterfall.
2. Stroke the strands to demonstrate that they are safe to touch.
3. Fan out the bright tips to cast a coloured light across a hand or face.
4. Spread the strands over a mat or seating and lay or sit the user on top.
5. Wrap the spray round limbs or the whole body.
6. Weave the strands in and out of a jumper or weave them in and out of netting which can then be draped over the user or suspended from a wall.
7. Use interactively with a switching system (ROMPA® Wi Fi or Sensory Magic™ for example).

## **11. Safety**

- Supervise the use of this product.
- Fibre Optic products are a constant source of safe visual and tactile stimulation but they should always be used under professional supervision. Fibre Optics are safe when correctly used and supervised. ROMPA cannot accept any liability for injury or damage caused by the abuse or misuse of these products. Though safe to stroke across the face, the strands should never be put in the mouth. These fibre optics are extremely tough and resistant to breaking if bitten, but nonetheless users should not be encouraged to mouth them.

## **12. Technical Specification**

Boss: 30mm in diameter  
Fibres: These plastic fibre optics are slightly thicker than standard, glass filled fibre optics  
Weight: e.g. 17903 – 1kg; 17905 – 2kg; 17907 – 3kg; 17908 – 4kg

## **13. Installation**

1. Remove the packaging and temporary ties that are holding the strands together.
2. Attach the boss to the light source – insert the boss (metal cylinder opposite end to the tips) into the port (cylindrical metal tube) of the light source.
3. Please see the separate instructions for your light source.

## **14. Care and Maintenance**

- Periodically wipe the fibre optics with a clean, slightly damp cloth. Switch off the light source first. Use a mild detergent and a slightly damp cloth, taking care not to wet the boss. Dry thoroughly before re-use. Do not immerse the fibre optics in water.
- Do not walk on the fibre optics or allow them to be repeatedly squashed, pinched, completely bent etc – once flattened or completely bent, the light will stop at that point and is unlikely to pass to the tips of the fibre optic strand.

## 15. Troubleshooting

*Please refer to the separate instructions for the Light Source*

Further copies of these instructions can be downloaded at [www.rompa.com](http://www.rompa.com)

We hope you find our instructions invaluable. If you have any suggestions for improving them further your comments will be greatly received – please contact us at [producthelp@rompa.com](mailto:producthelp@rompa.com)

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